

MULTIGATE Cadence® Hand Instrument



Clinch Grasper (Ratcheted)

Effortless precision

The Multigate Cadence® Clinch Grasper delivers effortless management of the operating field through precision and ergonomic design.

FEATURES

Comfort and Precision

Comfort and control are essential ingredients to a surgical procedure. Multigate's Cadence® Clinch Grasper delivers the ergonomic comfort and precision you need to perform at your best.

- ✓ Precision control
- ✓ Ergonomic handle design
- ✓ Ratchet on-off button
- ✓ Durable wide-opening jaws
- ✓ Recessed plug for improved safety
- ✓ Single-use eliminates CSSD costs



Ratcheted: Handle design incorporates the ability to grasp and lock tissue as desired using the on-off ratchet button.



Reduced fatigue: Ergonomic handle design provides optimal control during surgical procedure. Accommodates large and small hands, reducing fatigue.



Precision: Durable wide-opening jaws providing optimal application during the surgical procedure.

Cadence® Hand Instruments Clinch Grasper (Ratcheted)

RE-ORDER	SHAFT DIAMETER	W/L	QUANTITY
18-021L	Ø5mm	33cm	10pcs / box, 9boxes / carton
18-028L	Ø5mm	45cm	10pcs / box, 9boxes / carton

Designed for use with standard Australian reusable or single-use monopolar cables. Our Cadence® hand instruments have a monopolar 4mm banana-style plug which can accommodate cables with an outer plug diameter of up to 12mm. See IFU for more details.

For more information or to book a trial, contact your Multigate Territory Manager:
(02) 8718 2888 or visit [multigate.com.au/cadence](https://www.multigate.com.au/cadence)