

# MULTIGATE

## Cadence®

### Hand Instrument



## Johans Fenestrated Grasper (Ratcheted)

### Effortless precision

The Multigate Cadence® Johans Fenestrated Grasper delivers effortless management of the operating field through precision and ergonomic design.

#### FEATURES

### Comfort and Precision

Comfort and control are essential ingredients to a surgical procedure. Multigate's Cadence® Johans Fenestrated Grasper delivers the ergonomic comfort and precision you need to perform at your best.

- ✓ Precision control
- ✓ Ergonomic handle design
- ✓ Ratchet on-off button
- ✓ Durable wide-opening jaws
- ✓ Recessed plug for improved safety
- ✓ Single-use eliminates CSSD costs



**Ratcheted:** Handle design incorporates the ability to grasp and lock tissue as desired using the on-off ratchet button.



**Reduced fatigue:** Ergonomic handle design provides optimal control during surgical procedure. Accommodates large and small hands, reducing fatigue.



**Precision:** Durable wide-opening jaws providing optimal application during the surgical procedure.

#### Cadence® Hand Instruments Johans Fenestrated Grasper (Ratcheted)

RE-ORDER	SHAFT DIAMETER	W/L	QUANTITY
18-022L	Ø5mm	33cm	10pcs / box, 9boxes / carton
18-029L	Ø5mm	45cm	10pcs / box, 9boxes / carton

Designed for use with standard Australian reusable or single-use monopolar cables. Our Cadence® hand instruments have a monopolar 4mm banana-style plug which can accommodate cables with an outer plug diameter of up to 12mm. See IFU for more details.

For more information or to book a trial, contact your Multigate Territory Manager:  
**(02) 8718 2888** or visit [multigate.com.au/cadence](http://multigate.com.au/cadence)