

Trocar

Cannula Sets



Simple, safe and secure

The Multigate range of latex free trocar cannula sets provide simple, safe and secure access points for laparoscopic instruments to enter the abdominal cavity.

PRODUCT RANGE

- Autoshield Crown
- Bladeless Crown*
- Bladeless Optical Crown
- Secondary Bladeless
- Hasson Blunt Tip

FEATURES & BENEFITS

- ✓ Comprehensive range of styles and sizes including paediatric; 3 to 15 mm in diameter, 65 to 150mm in length
- ✓ Retractable valve design provides excellent sealing with low air leakage
- ✓ Detachable seal caps[^] for convenience of specimen extraction
- ✓ Sharp thin edged retractable blade[^] provides smoother puncture and helps prevent injury to organs
- ✓ Cannula shaft thread pattern[^] provides reliable abdominal wall fixation force to prevent accidental displacement of the cannula during surgery

[^] Not available on all trocar models



* Paediatric size available

 **MULTIGATE**

For more information or to book a trial, contact your Multigate Territory Manager:
(02) 8718 2888 or visit multigate.com.au

Trocar Cannula Sets



Simple, safe and secure

Trocar set with 1 x Trocar & 1 x Cannula (STERILE)

RE-ORDER	DESCRIPTION	SIZE (MM)	PIECES PER INNER	INNER PER CARTON
18-039	Bladeless Crown Trocar Set	3.0x65mm	6pcs / Inner box	18 boxes / ctn
18-040	Bladeless Crown Trocar Set	5.5x70mm	6pcs / Inner box	18 boxes / ctn
18-041	Bladeless Crown Trocar Set	5.5x100mm	6pcs / Inner box	18 boxes / ctn
18-042	Bladeless Crown Trocar Set	5.5x150mm	6pcs / Inner box	12 boxes / ctn
18-043	Bladeless Crown Trocar Set	10.5x110mm	6pcs / Inner box	18 boxes / ctn
18-044	Bladeless Crown Trocar Set	10.5x150mm	6pcs / Inner box	12 boxes / ctn
18-045	Bladeless Crown Trocar Set	12.5x110mm	6pcs / Inner box	18 boxes / ctn
18-046	Bladeless Crown Trocar Set	12.5x150mm	6pcs / Inner box	12 boxes / ctn
18-050	Bladeless Optical Crown Trocar Set	5.5x100mm	6pcs / Inner box	18 boxes / ctn
18-051	Bladeless Optical Crown Trocar Set	10.5x110mm	6pcs / Inner box	18 boxes / ctn
18-052	Bladeless Optical Crown Trocar Set	12.5x110mm	6pcs / Inner box	18 boxes / ctn
18-053	Bladeless Optical Crown Trocar Set	15x100mm	12pcs / Inner box	12 boxes / ctn
18-055	Autoshield Crown Trocar Set	5.5x100mm	6pcs / Inner box	18 boxes / ctn
18-056	Autoshield Crown Trocar Set	10.5x110mm	6pcs / Inner box	18 boxes / ctn
18-057	Autoshield Crown Trocar Set	12.5x110mm	6pcs / Inner box	18 boxes / ctn
18-060	Hasson Blunt Tip Trocar Set	12.5x110mm	6pcs / Inner box	18 boxes / ctn
18-061	Secondary Bladeless Trocar Set	5.5x100mm	15pcs / Inner box	18 boxes / ctn
18-063	Secondary Bladeless Trocar Set	5.5x150mm	12pcs / Inner box	12 boxes / ctn

Trocar set with 1 x Trocar & 2 x Cannula (STERILE)

RE-ORDER	DESCRIPTION	SIZE (MM)	PIECES PER INNER	INNER PER CARTON
18-047	Bladeless Crown Set with 1 x Trocar & 2 x Cannula	5.5x100mm	6pcs / Inner box	12 boxes / ctn
18-048	Bladeless Crown Set with 1 x Trocar & 2 x Cannula	10.5x110mm	6pcs / Inner box	12 boxes / ctn
18-049	Bladeless Crown Set with 1 x Trocar & 2 x Cannula	12.5x110mm	6pcs / Inner box	12 boxes / ctn
18-054	Bladeless Optical Crown Set with 1x Trocar & 2 x Cannula	5.5x100mm	6pcs / Inner box	18 boxes / ctn
18-058	Autoshield Crown Trocar Set with 1 x Trocar & 2 x Cannula	5.5x100mm	6pcs / Inner box	12 boxes / ctn
18-059	Autoshield Crown Trocar Set with 1 x Trocar & 2 x Cannula	12.5x110mm	6pcs / Inner box	12 boxes / ctn
18-062	Secondary Bladeless Trocar Set with 1x Trocar & 2 x Cannula	5.5x100mm	15pcs / Inner box	18 boxes / ctn