

MULTIGATE V-Text[®]

SURGICAL GOWNS

Uncompromising Safety and Performance.

V-Text[®] surgical gowns provide uncompromising performance and the reassurance of compliance with AAMI level 4 standards, making them perfect for long duration, high fluid level procedures.

TECHNOLOGY

Breathable viral barrier.

V-Text[®] gowns are made with a unique lightweight tri-laminate material that delivers a high performance impermeable but breathable viral barrier that allows the movement of warm air vapour away from the body whilst stopping bacteria and virus penetration.

FEATURES

AAMI level 4 protection.

Heat sealed seams and a responsive middle layer ensure our V-Text[®] gowns meet and exceed AAMI level 4 standards making them the ideal solution for high fluid, long duration procedures.

Comfort.

With generous dimensions to support upper body movement and a lightweight back panel to assist in regulating the wearer's body temperature, V-Text[®] gowns are designed to maximise comfort during procedures.

- ✓ Impervious
- ✓ AAMI level 4 protection
- ✓ Breathable
- ✓ Lightweight back panel
- ✓ Viral barrier
- ✓ Low linting

Standard Gown with 1 Huck Towel			
RE-ORDER	SIZE	STERILE / NON STERILE	CARTON QUANTITY
27-301	Medium	Sterile	20 Packs
27-305	Large	Sterile	20 Packs
27-306	X-Large	Sterile	20 Packs
27-307	XX-Large	Sterile	20 Packs
27-302	Large, X-Long	Sterile	20 Packs
27-303	X-Large, X-Long	Sterile	20 Packs
27-304	XX-Large, X-Long	Sterile	20 Packs



Lightweight
and breathable



Protection
against fluid
striekethrough
and viral
penetration

For more information please contact your Business Development Manager:
(02) 8718 2888 or visit multigate.com.au